Thank you for choosing VICTOR™, the all-round smart cane with LTE-M/NB-IoT technology. We appreciate your trust. Please reach out to us at victor@amptek-tech.com if you have any questions. Please visit our website at www.amptek-tech.com for news regarding our products.

VICTOR and the VICTOR logo are trademarks of Amptek Technologies Inc., registered in Canada and other countries. All rights reserved.

Specifications:

Model	VIC100
Height	76cm – 94cm
Weight	470g
Battery	3.70V 3000mAh Li-ion
LTE band	Cat-M1: B1, B2, B3, B4, B5, B8, B12, B13, B14, B17, B18, B19, B20, B25, B26, B28, B66 Cat-NB1/NB2: B1, B2, B3, B4, B5, B8, B12, B13, B17, B19, B20, B25, B26, B28, B66
GPS	L1 C/A supported
SIM	Nano SIM (included)
Charger	DC 5.0V 1000mA USB Type-A (included)
Weight capacity	124kg *
Handle material	PC-ABS & TPU / ABS & TPU
Tube material	6016-T6 aluminum alloy

* Support up to 25% of body weight.

HVIN: VIC100

PMN: VIC100

Contains FCC ID: 2BK6VVIC100

IC: 8115A-VIC100

Minimum radio transmission cycle occurs when a location request is received in which case approximately 7 kilobytes of data will be sent once every 5 minutes. This translates to approximately 0.056 seconds with LTE-M or 0.85 seconds with NB-IoT.

Features:



- Worldwide LTE-M/NB-IoT support with 1NCE nano SIM. Please visit www.1nce.com for a list of all supported countries and regions.
- GPS location with cellular location fallback.
- Panic button.
- Alert buzzer.
- Ambient light sensor.

- LED flashlight.
- Red/green color status indicator.
- Fall detection.
- Pedometer.
- Airplane mode for aviation safety.
- 3.7V, 3000mAh Li-ion battery with temperature sensor.
- USB Type-A magnetic charging cable and charger included.
- Sturdy 26mm diameter aluminum alloy tube.
- Replaceable rubber foot.
- Splash proof design.
- 8 adjustable heights ranging from 76cm to 94cm.
- Free automatic over-the-air firmware update.
- Free mobile app for 24/7 caring.
- First year wireless subscription included. Annual subscription renewal upon expiry.

Warning:


This product is intended to be used as a walking cane. This is not a toy. The product may only be used according to the instructions as described in this manual. All uses other than that described in this manual, are seen as unintended use. Do not attempt to disassemble the product. Modifications in the course of technological progress are reserved. We cannot be held in any case responsible for accidents or damage arising from the unintended use of the product. Moreover, any unintended use also causes the termination of the warranty period.

Walking pattern varies from time to time, person to person. VICTOR has the ability to continuously monitor and adapt to your pace. However, false alarm can happen. Each time a motion satisfies certain conditions to be considered as a fall, a doorbell sound can be heard. If VICTOR does not return to a straight upright position within 15 seconds, an alarm will be issued. Therefore, if you unintentionally triggered a doorbell sound, simply place VICTOR in a straight upright position for a second to cancel it.


Operations:

- Flip open the SIM cover and insert the supplied card tray and SIM card with contact pins facing down.

- VICTOR will start automatically once a valid SIM card is detected. It can take anywhere from 5 to 8 minutes for VICTOR to connect and register itself with the cloud server for the first time. The flashlight will remain on during the entire registration process. DO NOT remove the SIM card during the registration process.
- The flashlight will flash twice once VICTOR is connected with a wireless carrier and successfully registered with the cloud server.
- Download the mobile app from Google Play or App Store. Open the app and create an account. Enter the ICCID number labeled on the SIM card to connect with VICTOR. Please see **Mobile App** section for more details.
- Flashlight will turn on automatically when the ambient light is low and motion is detected. Flashlight will turn off automatically 10 seconds after motion stopped.
- Status light will flash red and green alternatively when battery voltage is low. Flip open the charger cover. Charge VICTOR using the supplied USB magnetic charging cable and the USB charger. Status light will remain in red during charging and will turn green once the battery is fully charged. All radio activities, including GPS, will be suspended during charging but will resume automatically once the charging cable is removed. One fully charged battery can operate approximately 7 days of normal usage.

- For safety reason, battery charging will be terminated automatically if the battery temperature rises to above 50°C.
- Press the panic button momentary for less than 5 seconds 3 times to enter airplane mode. Press the panic button momentary for less than 5 seconds 3 times again to exit airplane mode. The green status light flashes each time the panic button is being pressed. A double beep sound indicates entering airplane mode. The green status light will also flash for 5 seconds. A triple beep sound indicates exiting airplane mode. The flashlight will also flash for 5 seconds.
- Press and hold the panic button for 5 seconds or longer to issue an emergency alert. Alert will be sent to the mobile phone immediately. In addition, flashlight will flash and


buzzer will sound for help. Press and hold the panic button for 5 seconds or longer again to cancel the emergency alert.

- In case of a fall, an emergency alert will be sent to the mobile phone if the cane cannot be returned to a straight upright position within 15 seconds. Once an alert is issued, it can be canceled simply by returning the cane to a straight upright position.
- A warning message will be issued on the mobile app if no movement is detected from the cane for 8 consecutive hours. Please check the cane user immediately.
- Mobile app can monitor multiple canes simultaneously. Simply add the ICCID number of the SIM being used in the cane.
- VICTOR can be deactivated and removed from the cloud server through the mobile app. The red status light will flash once VICTOR is successfully deactivated and removed from the system. It is safe to remove the SIM card after VICTOR is decommissioned.

Mobile App:

Scan the QR code below to download the My Victor app from either the Apple App Store or the Google Play Store.

- Open the app and create an account.
- Tap on the blue [+] icon near the top right corner and enter the ICCID to pair with a cane. ICCID is the identification number of the SIM card used by the cane. It is a 19-digit number beginning with “89” printed on the SIM card. You can also tap on the barcode scanner icon  and then scan the barcode above the ICCID. Tap on the [+] icon again to pair with additional canes.

- Slide the user block to the right to edit the name or add a profile picture to identify the cane user.
- Slide the user block to the left to unpair or deactivate a cane. "Unpair" will only remove the cane from the mobile app. "Deactivate" will further delete the cane from the cloud server.
- Tap on [Get location] to see the location of the cane on a map. The location is updated once every 5 minutes when the map is displayed on the screen. A blue circle is used to indicate the accuracy of the current location. GPS location is only available when the cane is either in an outdoor environment or in indoor but close by a window. When GPS location is not available, the cane will fall back to use the cellular network for a coarse location.
- Pedometer shows the accumulated number of steps performed during a 24 hours period. The counter will be reset to zero everyday at midnight local time.

Warranty:

Each VICTOR carries a 1-year limited warranty. The warranty period commences upon the date of purchase by the original purchaser and continues for 12 months after that date. Warranty periods and warranty services are intended as a guide only and variations may occur.

The warranty covers manufacturing defects only. Please note that this does not include consumable items such as battery and rubber tip. The liability of Amptek Technologies is limited to the cost of repair and/or replacement of the product under warranty. The warranty is invalidated if the defect is caused (howsoever) by misuse, neglect, and tampering or incorrect adjustment. It is invalidated if unauthorized persons carry out any alterations and/or repairs. Also, the warranty is invalidated for repair due to not using the supplied charging accessories or nano SIM card.

To obtain technical assistance or to book a service/repair to your VICTOR under warranty, please contact us at **victor@amptek-tech.com**. To assist our customer service team, you will be asked to provide a copy of the original receipt as proof of purchase.

Product will be considered Out-of-Warranty and, if the product is deemed repairable, will be repaired at customer's sole cost and expense. Customer will be charged for all shipping and repair costs. Repair time can also be longer than usual. Detailed cost and estimated repair time can be confirmed by email.

All products will require verification and a determination whether it can be repaired. We reserve the right to refuse to repair the product at their discretion for any reason whatsoever.

Repair service may be refused due to the following reasons, including but not limited to:

- Non-genuine Amptek products.
- Customer misuse or abuse.
- Lack of replacement parts, tools or know-how.

Amptek Technologies does not guarantee that following the repair, the product will be fully functional according to the original specifications, liability or damages associated with such repairs.

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